GB

INSTRUCTIONS FOR USE

Silicone Foam Bordered Dressing

Farla

are

FarlaFOAM.S™

CE

0123

Description: Bordered Foam dressing with Silicone Gel Adhesive The dressing is made up of polyurethane film, hydrophilic polyurethane foam backing, super absorbent polyacrylate fibres, adhesive perforated silicone gel wound contact layer, and removable liners.

Indications: FarlaFOAM S is indicated for a wide range of exuding acute and chronic wounds including but not limited to pressure ulcers, leg/ foot ulcers and traumatic wounds e.g. skin tears and surgical wounds. The product can also be used on dry/necrotic wounds in combination with gels. FarlaFOAM S may be used as part of a prophylactic therapy to help prevent skin damage e.g. pressure, postoperative blistering.

Directions:

1. Cleanse the wound in accordance with normal procedures.
2. Dry the surrounding skin thoroughly.
3. Remove the removable backing liner and apply the adherent side to the wound. Do not stretch.
4. For best results, product should overlap the dry surrounding skin by at least 1-2 cm for the smaller sizes and 5cm for the larger sizes in order to protect the surrounding skin from maceration and fixate the dressing securely.
5. Product may be left in place for several days depending on the condition of the wound and surrounding skin, or as indicated by accepted clinical practice. Change at least every 7 days or according to the treatment plan.
6. Product can be used under compression bandaging.
7. Product can be used in combination with other dressings (see cautions) e.g. gels and gelling fiber dressings.

Contraindications: FarlaFOAM S should not be used on persons with a known allergy or sensitivity to the dressing or its components.

Cautions:

1. FarlaFOAM S should be disposed of in accordance with local procedures for infection control.
2. Do not re-use- this is a single use disposable dressing.
3. In the case of signs of clinical infection, consult a health care professional.
4. Do not use product together with oxidizing agents such as hypochlorite solutions or hydrogen peroxide.
5. The use of dressings as part of a prophylactic therapy does not preclude the need to continue to follow a comprehensive pressure ulcer prevention protocol, i.e. support surfaces, positioning, nutrition, hydration, skin care and mobility.

Sterilization: Sterilized using ethylene oxide.

Shelf Life: 3 years from date of manufacture, please see the product lot number and labels on the package.

Storage: This product should be kept dry and away from direct sunlight The dressing should be stored in dry conditions below 25°C

Symbols

	Batch code		Does not contain natural rubber latex
	Consult instructions for use		Do not re-use
	Do not use if package is damaged		Caution
	Sterilized using ethylene oxide		Upper lower limit of temperature
	Date of manufacture		Notified body number
	Use-by date		Keep Dry
	Keep away from sunlight		Manufacturer
			European Representative

Manufacturer: Huizhou Foryou Medical Devices Co., Ltd. North Shangxia Rd. Dongjiang Hi-tech Industry Park, 516005, Huizhou, P.R. China
www.foryoumedical.com

European representative: Shanghai International Holding Corp. GmbH (Europe), Eiffestrasse 80, 20537 Hamburg, Germany

Manufactured for and distributed by:
Farla Medical International Ltd. Staples Corner Business Park, Unit 2, 1000 North Circular Road LONDON NW2 7JP UK Tel: +44 (0)345 193 5 193
www.farlamedical.co.uk
www.farlainternational.com

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Rappresentante europeo: Shanghai International Holding Corp. GmbH (Europe), Eiffestrasse 80, 20537 Amburgo, Germania