



STANDING POSITION

[VIDEO](#)



Start by pulling the iD Care Net pants to mid-thigh and then turn the waistband down towards the knees to prepare to secure the pad. Fold the iD Form pad in half lengthways and channel it to activate the anti-leak cuffs.



Then pass the pad between legs from front to back with the back sheet turned away from the skin. Centre the pad by using the wetness indicators as a guide. Fan out both the back and the front of pad horizontally to create a channel shape between the legs. Then smooth out the product, front and back.



Finally, pull up the iD Care Net Pants into position.



Adjust accordingly.



LYING POSITION


[VIDEO](#)



Start by pulling the iD Care Net Pants to mid-thigh, then turn the waistband down towards the knees to prepare to secure the pad. Fold the iD Form pad in half lengthways to activate the anti-leak cuffs.



Gently roll the patient onto their side and, with safety in mind, bend the knees to create a channel between the legs. Then pass the pad between legs from front to back with the back sheet turned away from the skin and place the smallest area of the pad at the front. Centre the pad by using the wetness indicators as a guide. Fan out front and back making sure the anti leak cuffs stay activated.



Finally, pull up the iD Care Net Pants into position.



Adjust accordingly.