EN	A2 Instructions For Use Single Use Rongeurs Single Use Sterile Device, sterilised using Ethylene Oxide	EN
AR	A4 تعليمات الاستخدام قراضة العظم للاستخدام الواحد أداة معقمة للاستخدام لمرة واحدة، معقمة باستخدام أكسيد الإيثيلين	AR
ZH	A6 使用說明單次使用咬骨鉗 單次使用無菌器材，使用環氧乙烷消毒	ZH
DA	A8 Brugsanvisning Engangsrongeur Sterilt engangsudstyr, steriliseret med ethylenoxid	DA
NL	A10 Gebruiksaanwijzing Stans voor eenmalig gebruik Steriel hulpmiddel voor eenmalig gebruik, gesteriliseerd met ethyleenoxide	NL
FI	A12 Käyttöohjeet Kertäkäyttöinen luunpurija Kertäkäyttöinen steriili laite, steriloitu eteenioksidilla	FI
FR	A14 Mode d'emploi Pincés gouges à usage unique Dispositif stérile à usage unique, stérilisé à l'oxyde d'éthylène	FR
DE	A16 Gebrauchsanweisung Knochenstanzen zum Einmalgebrauch Steriles Einmalprodukt; mit Ethylenoxid sterilisiert	DE
EL	A18 Οδηγίες χρήσης Λαβίδες οστεοτομίας μίας χρήσης Στείρο προϊόν μίας χρήσης, αποστειρωμένο με οξειδίο του αιθυλενίου	EL
IS	A20 Notkunarleiðbeiningar Einnota beinaklippur Einnota sæft tæki, sæft með etýlenoxíði	IS
IT	A22 Istruzioni per l'uso Pinze ossivore monouso Dispositivo sterile monouso, sterilizzato con ossido di etilene	IT
NO	A24 Bruksanvisning Stansetang til engangsbruk Sterilt instrument til engangsbruk, sterilisert med etylenoksid	NO
PL	A26 Instrukcja użycia Odgryzacz do jednorazowego użytku Sterylny wyrób jednorazowego użytku, sterylizowany za pomocą tlenu etylenu	PL
PT	A28 Instruções de utilização Pinça goiva de utilização única Dispositivo estéril de utilização única, esterilizado por óxido de etileno	PT
SL	A30 Navodila za uporabo Klešče tipa rongeur za enkratno uporabo Sterilni pripomoček za enkratno uporabo, steriliziran z etilenoksidom	SL
ES	A32 Instrucciones de uso Pinzas gubias de un solo uso Producto estéril de un solo uso, esterilizado con óxido de etileno	ES
SV	A34 Bruksanvisning Engångsrongeur Steril engångsenhet, steriliserad med etylenoxid	SV

Products covered:

Kerrison Rongeurs - S125 Kerrison Rongeur 1mm x 200mm 45°; S126 Kerrison Rongeur 2mm x 200mm 45°; S127 Kerrison Rongeur 3mm x 200mm 45°; S137 Kerrison Rongeur 4mm x 200mm 45°; S138 Kerrison Rongeur 5mm x 200mm 45°; S161 Kerrison Rongeur 1mm x 200mm 45° for Delicate Use Only; S261 Kerrison Rongeur 2mm x 200mm 45° for Delicate Use Only; S262 Kerrison Rongeur 3mm x 200mm 45° for Delicate Use Only; S367 Kerrison Rongeur 2mm x 170mm 90°.

Intended Purpose:

Instrument used for resection of bone.

Indication for Use:

Clinical Indications for the Surgical Kerrison devices are for procedures where removal of bone is required.

Intended Clinical Benefits:

Facilitates precision removal of bone fragments to allow surgical access and/or alleviate pressure on surrounding tissue.

Providing the user with the option of a single use instrument to reduce the likelihood of cross contamination and procedural site infection when compared to reusable alternatives.

Intended Patient Groups:

Surgical Kerrison Devices are intended for use in general adult and paediatric patient populations where resection of bone is required. Patient suitability is at the discretion of a trained healthcare professional.

Contraindications:

Not for use in **direct contact** with the central nervous system as defined in EU MDR 2017/745.

Not for use in **direct contact** with the central circulatory system as defined in EU MDR 2017/745.

Intended User Qualification and Environment of Use:

For use by trained healthcare professionals in a surgical setting.

⚠Warnings:

- **Do not** use expired product.
- **Do not** use if device or sterile barrier is damaged.
- **Do not** apply excessive force to the device as this could cause device to break or injure patient/user.
- **Do not** attempt to resect large sections of bone with one activation of the device. Doing so may result in reduced device performance or incomplete resection of target bone, reduced precision of resection, patient harm and/or damage to the device.
- **Do not** use in conjunction with MRI (magnetic resonance imaging).
- Use in line with sterile practice guidelines.
- Device is for single use only and cannot be reused or reprocessed due to risks associated with infection and cross contamination, potential damage to the device and patient injury.
- Use the device according to its intended purpose. Use of the device to resect soft tissue may result in reduced performance.
- Always select the appropriate device size. If a device is selected with a footplate too large for the target resection site, this may result in tissue damage.
- Avoid overstraining the product. If target bone is found to not be cleanly resected after activation of the device, avoid twisting, pulling or levering the device to complete resection. Doing so may result in reduced precision of resection, patient harm and/or damage to the device.
- Bone fragment **MUST** be removed from the jaw of device after every cut. Failure to do so may result in reduced performance of the device.
- Use with care in conjunction with heat sources or lasers to avoid injury to surgical staff or patient.

Precautions:

- Inspect all sterile packaging prior to removing the device from packaging and dispose of device if sterile barrier is broken or damaged.
- Inspect device for damage prior to use.
- Confirm suitability of device geometry prior to use. To aid visual identification, the bite size of the device is etched on the side of Surgical Kerrison Devices.
- Use with care in conjunction with devices such as electrosurgical devices, heat sources, power tools, or lasers to avoid injury to user or patient, or damage to the device.
- If device becomes damaged during use, discontinue to use the device.
- 'Delicate Use Only' (DUO) variants are available with low profile footplates for use in small spaces.
- The jaw of the device contains sharp edges. Take care while handling to avoid injuries to users or patients.
- Device mechanism contains moving parts. Take care while handling and operating to avoid trapping hand/glove.

Use Instructions:

1. Read and understand the instruction for use before using the device.
2. Confirm device specification and geometry is suitable for clinical purpose.
3. Inspect all sterile packaging prior to removing the device from packaging and dispose of device if sterile barrier is broken or damaged.
4. Inspect label for expiry date and do not use expired product.
5. Remove and handle the device in accordance with sterile practice.
6. Inspect device for damage prior to use, damaged devices should not be used. Unused, damaged items should be reported and returned to the manufacturer.
7. Perform functional check prior to use.
Introduce device to procedure site and navigate jaw aperture to target bone under vision, ensuring footplate of device does not apply excessive force to surrounding tissues.
8. When satisfied the jaw aperture is safely positioned, perform a single cut.
9. After each cut ensure bone fragment is removed from the jaw of the device.
10. Check device is intact following completion of the procedure.

Storage and Disposal:

- Store in a dry environment within temperature limits and shelf life as marked. Do not open until required.
- Dispose of device as clinical waste or sharps in accordance with clinical waste guidelines.

Manufacturing Details:

Contact the manufacturer for technical support. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority.

Symbol Key:

	Manufacturer		Caution
	Catalogue number		Medical device
	Use-by-date		Unique device identifier
	Date of manufacture (YYYY-MM-DD)		Do not resterilize
	Authorised representative in the european community/european union		CAUTION: US Federal law restricts this device to sale by or on the order of the physician
	Batch code		Quantity of content
	Single sterile barrier system with protective packaging outside		Indicates recovery/recyclable package or package material
	Do not reuse		Latex Free
	Do not use if packaging is damaged and consult instructions for use		Consult instructions for use
	Sterilised using ethylene oxide		Single sterile barrier system
	Importer		Temperature limit
	Keep dry		UKCA certification mark
	Authorised representative in Switzerland		CE certification mark

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Advena Ltd, Tower Business Centre, 2nd Flr.,