

# clinell®

## Alcohol Wipes

---

Surface disinfection  
of non-invasive  
medical devices.



# Rapid, high quality alcohol wipes

## Ultra-quick drying, pharmaceutical grade 70% isopropyl alcohol wipes.

Disinfection wipes, which leave no smear, for surface disinfection of non-invasive medical devices.

### Proven efficacy



High-quality durable wipes containing pharmaceutical grade 70% isopropyl alcohol. Our wipes are proven to have bactericidal and virucidal action.

### Unique packaging



Wipes are dispensed one at a time, so the top of the next wipe is not dry, which commonly occurs with tub dispensers. This eliminates the infection control risk of dry contaminated wipes.

### Less plastic. Less waste



Our flow wrap packaging is environmentally friendly as there is much less packaging to dispose of compared to standard alcohol tubs.

### Superior quality



Clinell Alcohol Wipes contain 25gsm non-woven spun bond, which is wetter & thicker than competing alcohol wipes allowing for larger area disinfection with fewer wipes, saving time and money.

### For the surface disinfection of non-invasive medical devices:

- Ultrasound probes
- Pulse oximeter probes
- Doppler probes
- Stethoscopes
- ECG machines
- Blood pressure machines
- Blood pressure cuffs
- Radiology equipment
- Incubators
- Monitors

# Antimicrobial efficacy





Proven to kill at least 99.999% of microorganisms.

Including MRSA, *E. coli*, *Pseudomonas*, *Enterococcus*, *Candida* and *Aspergillus* according to EN1276, EN13727, EN14561, EN13624 & EN14562.

	Microorganism example	EN test
<b>Bacteria</b>	<i>Acinetobacter baumannii</i>	EN13727
	<i>Enterococcus hirae</i>	EN13727
		EN 16615
	<i>Escherichia coli</i> K12	EN13727
	Methicillin resistant <i>Staphylococcus aureus</i>	EN13727
	<i>Pseudomonas aeruginosa</i>	EN13727
		EN16615
	<i>Staphylococcus typhimurium</i>	EN13727
	<i>Staphylococcus aureus</i>	EN13727
		EN16615
<b>Fungi</b>	<i>Aspergillus brasiliensis</i>	EN13624
<b>Mycobacteria</b>	<i>Mycobacterium avium</i> (2 year aged)	EN14348
	<i>Mycobacterium terrae</i> (2 year aged)	EN14348
<b>Viruses</b>	Norovirus	EN14476
	Murine coronavirus	EN14476
	Vaccinia virus	EN14476
<b>Yeast</b>		EN13624
	<i>Candida albicans</i>	EN16615

# Ordering information

Alcohol Wipes - Pack of 150		Transport unit view
	Product code	CAW150L
	Wipes per pack	150
	NHSSC code (UK only)	VJT839
	Expiry date	2 years
	Unit per carton	18
	Carton size	390 x 330 x 190 mm
	Carton per layer	9
	Layer height	1284 mm
	Cartons per pallet	54
	Classification	Class IIa Medical Device CE Mark
		

For more information, speak to your GAMA Healthcare Area Manager or visit [www.clinell.com](http://www.clinell.com)

**CE 0050** Class IIa Medical Device  
Conforms to the Medical Device Directive, 93/42/EEC.

Use disinfectants safely. Always read the label and product information before use. Always follow the medical device manufacturer's cleaning and disinfection guidelines prior to use with Clinell Alcohol Wipes.



**GAMA Healthcare Ltd.,**  
The Maylands Building, Maylands Avenue,  
Hemel Hempstead, HP2 7TG, UK.

T: +44 (0)20 7993 0030  
E: [info@gamahealthcare.com](mailto:info@gamahealthcare.com)  
[www.gamahealthcare.com](http://www.gamahealthcare.com)

JBN220359