

Sterilisation Pouches and Sterilisation Rolls (52***)

Pre-formed sterile barrier systems for the safe packaging of medical devices and instruments prior to terminal sterilisation

Technical Data:

Sterilisation method: Steam (STEAM) and EO
 Paper weight: approx. 60 g/m²
 Sealing temperature: 150–180°C
 Shelf life: 5 years from date of manufacture

Standards and Compliance:

Complies with MDR 2017/745, EN ISO 11607-1/-2, EN 868-3, EN 868-5 and EN ISO 11140-1 Type 1. Quality management according to EN ISO 9001 and EN ISO 13485.

Storage:

Store in a cool, dry place. Protect from direct sunlight and heat sources.

For more information about the range click on the Image 



Product Code and Size:

Code	Size
52100/ RP050200	50 x 200 mm Sterilisation Roll
52102/ RP075200	75 x 200 mm Sterilisation Roll
52104/ RP100200	100 x 200 mm Sterilisation Roll
52106/ RP 150200	150 x 200 mm Sterilisation Roll
52108/ RP200200	200 x 200 mm Sterilisation Roll
52090/ BA090250	90 x 250 mm Sterilisation Pouch
52140/ BA140260	140 x 260 mm Sterilisation Pouch

Product Description

Our sterilisation pouches and sterilisation rolls are pre-formed sterile barrier systems for the safe packaging of medical devices and instruments prior to terminal sterilisation. The products are suitable for steam sterilisation (autoclave) and ethylene oxide (EO) sterilisation.

Features

- PET/PP laminated film
- Medical-grade paper (approx. 60 g/m²)
- Integrated sterilisation indicators
- Latex-free
- DEHP-free
- High tear resistance
- Optimum visibility of the contents


Benefits

- Suitable for steam and EO sterilisation
- Complete batch traceability
- High material strength
- Easy handling
- A wide range of sizes and formats

Applications

- Hospitals
- Clinics
- Medical and dental practices
- Laboratories
- Medical device manufacturers
- Industrial sterilisation processes

Sterilisation Roll



Sterilisation Pouch



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com

