



Contact® (150*)

offline washed · chemical tested · white · medical grade
medium weight · latex · examination gloves

Physical Specification:

Material: Latex (Flat)
Colour: White
Weight: 6.2 g
Thickness: 0.11 mm
Finish: Powder-free
Fit: Ambidextrous
Sterile: No

Order Codes | Sizing:

1501 | Extra Small
1502 | Small
1503 | Medium
1504 | Large

Packaging:

Each case contains 10 boxes
Each dispenser box contains 100 pairs

Classification:



For more about the **Premium** range, please **click on the image**



Features

- Washed offline
- Tested to EN374 Type B
- AQL 1.5
- Free from powder
- Tested to EN1186
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Resistant to a range of chemicals
- Meets medical standards for pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves
- Extremely grippy and reliable on the outside, the glove feels soft against the skin on the inside and slides effortlessly onto the hand

Sectors

- Surgery
- Laboratories
- Clinics
- Chemical manufacture and use
- Dentistry
- Dental hygiene
- Emergency services



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
de.unigloves.com

