

# Product specification sheet

Revision date 18th September 2025

|     |                          |                                           |         |                  |         |              |                          |           |         |         |               |
|-----|--------------------------|-------------------------------------------|---------|------------------|---------|--------------|--------------------------|-----------|---------|---------|---------------|
| 1.  | Product code             | 600X*                                     |         |                  |         |              |                          |           |         |         |               |
| 2.  | Product description      | Format Blue                               |         |                  |         |              |                          |           |         |         |               |
| 3.  | Material                 | Nitrile                                   |         |                  |         |              |                          |           |         |         |               |
| 4.  | Inspection parameters    | Criteria                                  |         | Inspection level |         | AQL          |                          | Standard  |         |         |               |
|     |                          | Dimensions                                |         | S-2              |         | 4.0          |                          | EN455-2   |         |         |               |
|     |                          | Physical properties                       |         | S-2              |         | 4.0          |                          | EN455-2   |         |         |               |
|     |                          | Water tight test (1,000ml)                |         | G-1              |         | 1.5          |                          | EN455-1   |         |         |               |
|     |                          | Visual defects                            |         | G-1              |         | 2.5          |                          | In-house  |         |         |               |
| 5.  | Glove dimensions         | Dimensions                                |         | Size (mm)        |         |              |                          |           |         |         |               |
|     |                          |                                           |         |                  | S       | M            | L                        | XL        |         |         |               |
|     |                          | Product code                              |         |                  | 6002    | 6003         | 6004                     | 6005      |         |         |               |
|     |                          | Palm width                                |         |                  | 85±5    | 95 ± 5       | 106 ± 5                  | ≥120 ± 10 |         |         |               |
|     |                          | Length                                    |         |                  | ≥ 240   | ≥ 240        | ≥ 240                    | ≥ 240     |         |         |               |
| 6.  | Glove thickness          | Point of thickness measurement            |         |                  |         |              | Nominal single wall (mm) |           |         |         |               |
|     |                          | Cuff (at beaded end)                      |         |                  |         |              | 0.07                     |           |         |         |               |
|     |                          | Palm (at centre of palm)                  |         |                  |         |              | 0.09                     |           |         |         |               |
|     |                          | Fingertip (at fingertip)                  |         |                  |         |              | 0.13                     |           |         |         |               |
| 7.  | Physical properties      | Parameters                                |         | Before aging     |         | After aging  |                          |           |         |         |               |
|     |                          | Median force at break (N)                 |         | Min 6.0          |         | Min 6.0      |                          |           |         |         |               |
|     |                          | Tensile strength (MPa)                    |         | Min 14.0         |         | Min 14.0     |                          |           |         |         |               |
|     |                          | Ultimate elongation (%)                   |         | Min 500          |         | Min 400      |                          |           |         |         |               |
| 8.  | Proteins                 | None                                      |         |                  |         |              |                          |           |         |         |               |
| 9.  | Powder residue           | Maximum 2.0 mg / glove                    |         |                  |         |              |                          |           |         |         |               |
| 10. | Cuff type                | Beaded                                    |         |                  |         |              |                          |           |         |         |               |
| 11. | Design                   | Ambidextrous                              |         |                  |         |              |                          |           |         |         |               |
| 12. | Surface                  | Microroughened                            |         |                  |         |              |                          |           |         |         |               |
| 13. | Colour                   | Blue                                      |         |                  |         |              |                          |           |         |         |               |
| 14. | Quality assurance        | Manufactured under ISO 9001 and ISO 13485 |         |                  |         |              |                          |           |         |         |               |
| 15. | Weight                   | 4.4g (size medium)                        |         |                  |         |              |                          |           |         |         |               |
| 16. | Shelf life               | 5 years                                   |         |                  |         |              |                          |           |         |         |               |
| 17. | Standards compliant with | EN455-1                                   | EN455-2 | EN455-3          | EN455-4 | EN ISO 21420 | EN374-1                  | EN374-2   | EN374-4 | EN374-5 | Food approved |
|     |                          | ✓                                         | ✓       | ✓                | ✓       | ✓            | ✓                        | ✓         | ✓       | ✓       | ✓             |

The above data represents the product norm, which may be changed without notice and is intended as a guide only.