



Format Blue (600X*)

offline washed · chemical tested · blue · medical grade
medium weight · nitrile · examination gloves

Physical Specifications:

Material: Nitrile (latex free)
Colour: Blue
Length: 240 mm
Weight: 4.4 g
Thickness: 0.09 mm
Finish: Powder free, chlorinated
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

6002 | Small
6003 | Medium
6004 | Large
6005 | Extra large

Packaging:

Each dispenser box contains 100 gloves
Each case contains 10 dispenser boxes

Classification:



For more about the **Premium** range, please click on the image



Features


- Washed offline
- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace
- Chemical manufacture and use
- Dentistry • Dental hygiene
- Electronics • Engineering
- Food processing
- Healthcare
- Laboratory work
- Manufacturing and assembly
- Pharmaceutical work
- Printing



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
de.unigloves.com



Issue No. 2

Prepared at Gillingham
on 29/09/25