

Product specification sheet

Revision dated May 2022

1.	Product code	GA006*									
2.	Product description	PRO-TECT Black HD									
3.	Material	Nitrile									
4.	Inspection parameters	Criteria		Inspection level	AQL	Standard					
		Dimensions		S-2	4.0	EN455-2					
		Physical properties		S-2	4.0	EN455-2					
		Water tight test (1,000ml)		G-1	1.5	EN455-1					
		Visual defects		G-1	2.5	In-house					
5.	Glove dimensions	Dimensions	Size (mm)								
			XS	S	M	L	XL				
		Product code	GA0061	GA0062	GA0063	GA0064	GA0065				
		Palm width	≤80	80 ± 10	95 ± 10	110 ± 10	≥ 110				
		Length	≥ 240	≥ 240	≥ 240	≥ 240	≥ 240				
6.	Glove thickness	Point of thickness measurement			Nominal single wall (mm)						
		Cuff (at beaded end)			0.09						
		Palm (at centre of palm)			0.12						
		Fingertip (at fingertip)			0.15						
7.	Physical properties	Parameters		Before aging	After aging						
		Median force at break (N)		Min 6.0	Min 6.0						
		Tensile strength (MPa)		Min 14.0	Min 14.0						
		Ultimate elongation (%)		Min 500	Min 400						
8.	Proteins	None									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Microroughened									
13.	Colour	Black									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	6.2g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.