

## Pearl (GP\*\*\*\*)

16 colours · medical grade · medium weight  
nitrile · examination gloves

### Physical Specifications:

Material: Nitrile (latex free)  
Colours: Black, Blue, Burgundy, Cobalt, Green, Lime, Magenta, Mint, Opal, Orange, Pink, Red, Sapphire, Violet, White, Yellow  
Length: 240 mm  
Weight: 3.5 g  
Thickness: 0.07 mm  
Finish: Powder free, chlorinated  
Fit: Ambidextrous  
Sterile: No

### Order Codes and Size Information:

GP\*\*\*1 | Extra small  
GP\*\*\*2 | Small  
GP\*\*\*3 | Medium  
GP\*\*\*4 | Large  
GP\*\*\*5 | Extra large

GP003X Black / GP001X Blue  
GP012X Burgundy / GP101X Cobalt  
GP004X Green / GP350X Lime  
GP850X Magenta / GP680X Mint  
GP009X Opal / GP013X Orange  
GP005X Pink / GP006X Red  
GP008X Sapphire / GP007X Violet  
GP002X White / GP011X Yellow

**Note:** Not all colours come in Extra Small

### Packaging:

Each dispenser box contains 100 gloves

Each case contains 10 dispenser boxes

### Classification:



For more about the Pearl range, please click on the image



### Features

- Tested to EN374
- Available in sixteen colours
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Touchscreen compatible

### Benefits

- Resistant to a range of chemicals
- Ideal for colour-coding in food, cleaning and healthcare environments
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

### Sectors

- Care home work
- Chemical manufacture and use
- Cleaning and janitorial
- Dentistry • Dental hygiene
- Emergency services
- Floristry and gardening
- Food processing • Hair and beauty
- Healthcare • Hospital and clinical duties
- Janitorial • Painting and decorating
- Piercing • Tattoo

### Colour range:



**UNIGLOVES UK**  
Registered in England No: 4010200  
VAT registration 791817787

3 Ambley Green, Gillingham Business Park,  
Gillingham, ME8 0NJ, UK  
+44 (0)800 049 6602  
enquiries@unigloves.co.uk  
uk.unigloves.com

**UNIGLOVES GERMANY**  
Arzt- & Klinikbedarf  
Handelsgesellschaft mbH

Camp-Spich-Straße 71  
D-53842 Troisdorf-Spich,  
+49 (0) 2241-9323-0  
unigloves@unigloves.de  
unigloves.de

