

# Product specification sheet

Revision dated May 2022

1.	Product code	GS009*									
2.	Product description	Unicare Blue Vinyl Powdered									
3.	Material	Vinyl									
4.	Inspection parameters	Criteria		Inspection level	AQL	Standard					
Dimensions		S-2	4.0	EN455-2							
Physical properties		S-2	4.0	EN455-2							
Water tight test (1,000ml)		G-1	4.0	EN455-1							
Visual defects		G-1	2.5	In-house							
5.	Glove dimensions	Dimensions		Size (mm)							
				XS	S	M	L	XL			
Product code		-	GS0092	GS0093	GS0094	GS0095					
Palm width		-	80 ± 10	95 ± 10	110 ± 10	≥ 110					
Length		-	≥ 240	≥ 240	≥ 240	≥ 240					
6.	Glove thickness	Point of thickness measurement				Nominal single wall (mm)					
Cuff (at beaded end)		0.05									
Palm (at centre of palm)		0.06									
Fingertip (at fingertip)		0.07									
7.	Physical properties	Parameters		Before aging		After aging					
Median force at break (N)		N/A		N/A							
Tensile strength (MPa)		Min 9		Min 9							
Ultimate elongation (%)		N/A		N/A							
8.	Proteins	None									
9.	Powder residue	Maximum 10.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Smooth									
13.	Colour	Blue									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	4.0g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
						✓				✓	

The above data represents the product norm, which may be changed without notice and is intended as a guide only.