

Product specification sheet

Revision dated July 2026

1.	Product code	GT001*									
2.	Product description	SELECT Black 300									
3.	Material	Latex									
4.	Inspection parameters	Criteria		Inspection level		AQL		Standard			
		Dimensions		S-2		4.0		EN455-2			
		Physical properties		S-2		4.0		EN455-2			
		Water tight test (1,000ml)		G-1		1.5		EN455-1			
		Visual defects		G-1		2.5		In-house			
5.	Glove dimensions	Dimensions		Size (mm)							
				XS		S	M		L	XL	
		Product code		GT0011	GT0012	GT0013		GT0014		GT0015	
		Palm width		≤80	80 ± 10	95 ± 10		110 ± 10		≥ 110	
		Length		≥ 300	≥ 300	≥ 300		≥ 300		≥ 300	
6.	Glove thickness	Point of thickness measurement					Nominal single wall (mm)				
		Cuff (at beaded end)					0.08				
		Palm (at centre of palm)					0.13				
		Fingertip (at fingertip)					0.16				
7.	Physical properties	Parameters		Before aging			After aging				
		Median force at break (N)		Min 6.0			Min 6.0				
		Tensile strength (MPa)		Min 18.0			Min 14.0				
		Ultimate elongation (%)		Min 650			Min 500				
8.	Proteins	< 50 µg/g extractable latex proteins									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Microroughened									
13.	Colour	Black									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	7.8g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN ISO 21420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.