

Format White (100*)

Disposable Examination Gloves made of Nitrile

Nitrile • Latex free • Powder free • White • Non-sterile •

Chemical protective glove EN 374 • Suitable for short-term food contact •

Protection against viruses, fungi and bacteria



Specifications:

Material: Nitrile
 Powder: Powder free
 Ingredients: Thiuram free, no proteins contained
 Cuff: Beaded
 Design: Ambidextrous
 Sterile: Non-sterile
 Design: Micro-roughened fingers and palm
 Colour: White
 Quality control :
 EN 374 CAT III - Typ C;
 EN ISO 21420;
 EN 455 Teil 1-4;
 DIN EN ISO 13485:2021

Item number and size:

1001 | Extra Small
 1002 | Small
 1003 | Medium
 1004 | Large
 1005 | Extra Large

Packaging:

100 pieces per box
 10 boxes per carton

Classification:



Test Parameter

Test Chemical	Inspection Level	AQL
Dimensions	EN 455-2	Median value meets the requirements
Physical Properties	EN 455-2	Median value meets the requirements
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	2.5
Minor Visual defects	G-1	4.0

Dimensions

Item no.	1001	1002	1003	1004	1005
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	75 ± 5 mm	85 ± 5 mm	95 ± 5 mm	106 ± 5 mm	120 ± 10 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer (mm)
Cuff (at a distance of 40mm to the beaded edge)	0.07 ± 0.02
Palm (in the middle of the palm)	0.09 ± 0.03
Fingertips (at a distance of 15mm to the tip of the middle finger)	0.12 ± 0.03

Physical Properties

Parameter	Before Aging	After Aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

CCQS Certification Services Limited, Block 1 Blanchardstown Corporate Park,
 Ballycoolin Road, Blanchardstown, D15 AKK1, Ireland
 (Notified Body 2834)

Test chemical	EN ISO 374-1:2016 + A1:2018 Permeation level	EN ISO 374-4:2019 Degradation %
Sodium hydroxide 40% (K)	6	-5.3

Hint

The above information apply to the standard version.
 Subject to change.



UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
 D-53842 Troisdorf
 ☎ +49 (0) 2241-9323-0
 ✉ unigloves@unigloves.de
 🌐 unigloves.de

