

Unigloves Standard Plus (105X*)

natural · medical grade · light weight
 latex · examination gloves

Physical Specification:

Material: Latex (Flat)
 Colour: White
 Thickness: 0.08 mm
 Finish: Powderfree
 Fit: Ambidextrous
 Sterile: No

Order Codes | Sizing:

1051 | Extra Small
 1052 | Small
 1053 | Medium
 1054 | Large
 1055 | Extra Large

Packaging:

Each case contains 10 boxes
 Each dispenser box
 contains 100 pairs

Classification:



For further
 information
 please click
 on the
 image



Features

- AQL 1.5
- Free from powder
- Free from phthalates
- Tested to EN1186
- Tested to EN374 Type B
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Meets medical standards for pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with food
- Suitable for contact with all food types
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Care home work
- Cleaning and janitorial
- Dentistry
- Dental hygiene
- Emergency services
- Food processing
- General assembly
- Healthcare
- Hospital and clinical duties
- Janitorial
- Local Authority work
- Painting and decorating
- Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com



Issue No. 1

Prepared at Gillingham
 on 18/11/2025