

Standard Plus (105*)

Disposable Examination gloves made of latex

Latex • Powder free • Natural colour • Non-sterile • Chlorinated •

Chemical protective glove EN 374 • suitable for short-term food contact •

Protection against viruses, fungi and bacteria

Specifications:

Material: High quality natural latex

Powder: Powder free

Protein content: max.
50 µg/dm²

Cuff: Beaded cuff

Design: Ambidextrous

Surface: Micro-roughened fingertips

Colour: Natural colour

Quality control:

EN 374 CAT III - Type B,
EN 21420;

EN 455 Part 1-4;

EN ISO 13485:2016

Item number and size:

1051 | Extra Small

1052 | Small

1053 | Medium

1054 | Large

1055 | Extra Large

Packaging:

100 pieces per box

10 boxes per carton

Classification:



Für weitere
Informationen
klicken Sie
bitte auf
das Bild



Test parameter

Test criteria	Level	AQL
Dimension	S-2	4,0
Physical properties	S-2	4,0
1000-ml-Watertightness test	G-1	1,5
Visual inspection for major errors	G-1	2,5
Visual inspection for minor errors	G-1	4,0

Dimensions

Item.-no.	1051	1052	1053	1054	1055
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	75 ± 5 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring points	One layer (mm)
Palm (in the middle of the palm)	min. 0,08 mm

Physical properties

Parameter	Before aging	After aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

SATRA Technology Europe, Bracetown Business Park, Clonee, D15 YN2P Ireland
(Notified Body 2777)

Test chemical	EN ISO 374-1:2016+A1:2018/Type B Permeation Level	EN ISO 374-4:2019 Degradation %
Sodium hydroxide 40% (K)	4	-72.4
Sulphuric acid 96% (L)	1	92.9
Nitric acid 65% (M)	1	20.9
Formaldehyde 37% (T)	6	-81.7
Hydrogen peroxide 30% (P)	6	-64.1

Hint

The above information applies to the standard version.
Subject to changes.



UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
D-53842 Troisdorf
☎ +49 (0) 2241-9323-0
✉ unigloves@unigloves.de
🌐 unigloves.de

SUBSIDIARY OF THE
UNIGLOVES®
UG HEALTHCARE GROUP
ughealthcarecorporation.com