

Unigloves Expert Plus (12XX*)

Surgical disposable gloves made of latex · powder-free
white · sterile · anatomical fit

Physical Specification:

Material: Latex (Flat)
 Colour: Natural
 Length: 295 mm
 Thickness: 0.22 mm
 Finish: Powder-free
 Fit: Ambidextrous
 Sterile: Yes

Order Codes | Sizing:

1260 | Extra Small
 1265 | Small
 1270 | Medium
 1275 | Large
 1280 | Extra Large
 1285 | Extra Extra Large

Packaging:

Each case contains 4 boxes
 Each dispenser box contains 50 pairs

Classification:



AQL 1.5 (Pinholes)

EN455 (Medical Grade)



For further
 information
 please click
 on the
 image



Features

- Sterile disposable glove made from high-quality natural latex
- Anatomical fit
- Micro-roughened palm for improved grip
- High tear resistance despite elasticity
- Low protein content ($<50 \mu\text{g}/\text{dm}^2$), skin-friendly
- Sizes: XS to L
- Powder-free, sterile packaging in pairs

Benefits

- Maximum safety and hygiene in the operating room
- Precise movements and superior comfort
- Secure handling of instruments, even in wet conditions
- Durable and reliable for all surgical procedures
- Gentle on sensitive skin
- Suitable for all medical and dental professionals
- No compromise on quality and comfort

Sectors

- Surgery
- Laboratories
- Clinics
- Chemical manufacture and use
- Dentistry
- Dental hygiene
- Emergency services



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com



Issue No. 1

Prepared at Gillingham
 on 18/11/2025