

# Expert Plus (12\*\*)

## Disposable Latex Examination Gloves

Latex • Powder free • White • Sterile • Anatomical fit

### Specifications:

Material: Natural High Grade Rubber Latex  
 Powder: Powder free (< 2 mg/glove)  
 Protein content: < 50 µg/g extractable latex proteins  
 Cuff: beaded cuff  
 Design: Hand-specific  
 Colour: White  
 Quality control:  
 EN 455 Teil 1-4;  
 DIN EN ISO 13485:2021

### Item Number and Size:

1260 | 6.0  
 1265 | 6.5  
 1270 | 7.0  
 1275 | 7.5  
 1280 | 8.0  
 1285 | 8.5

### Packaging:

50 pair per dispenser  
 4 dispensers per carton

### Classification:



### Test parameter

| Test criteria                      | Level    | AQL                                 |
|------------------------------------|----------|-------------------------------------|
| Dimensions                         | EN 455-2 | Median value meets the requirements |
| Physical Properties                | EN 455-2 | Median value meets the requirements |
| Freedom from Holes                 | G-1      | 1.5                                 |
| Visual inspection for major errors | G-1      | 1.5                                 |
| Visual inspection for minor errors | G-1      | 2.5                                 |

### Dimensions

| REF    | 1260      | 1265      | 1270      | 1275      | 1280       | 1285       |
|--------|-----------|-----------|-----------|-----------|------------|------------|
| Size   | 6.0       | 6.5       | 7.0       | 7.5       | 8.0        | 8.5        |
| Width  | 77 ± 5 mm | 83 ± 5 mm | 89 ± 5 mm | 95 ± 5 mm | 102 ± 6 mm | 108 ± 6 mm |
| Length | ≥ 295 mm  | ≥ 295 mm  | ≥ 295 mm  | ≥ 295 mm  | ≥ 295 mm   | ≥ 295 mm   |

### Material thickness

| Measuring points                                                   | Double Wall (mm) |
|--------------------------------------------------------------------|------------------|
| Cuff (at a distance of 25 mm to the beaded edge)                   | 0.18 ± 0.02      |
| Palm (in the middle of the palm)                                   | 0.22 ± 0.02      |
| Fingertips (at a distance of 15mm to the tip of the middle finger) | 0.23 ± 0.02      |

### Physical Properties

| Standard           | Before aging | After aging |
|--------------------|--------------|-------------|
| Force at Break (N) | min. 9.0     | min. 9.0    |

### Quality assurance

Comply the requirements of ASTM D3577-91, EN 455 part 1-4  
 Manufactured according GMP and EN ISO 9001 and EN 13485 quality assurance system

### Hint

The above information apply to the standard version.  
 Subject to change.



**TERANG NUSA Sdn. Bhd.**  
 1 Jalan 8  
 Pengkalan Chepa 2 Industrial Zone  
 16100, Kota Bharu, MALAYSIA

**EC REP**  
**ULMA INTERNATIONAL GmbH**  
 Pfaffenweg 35  
 89231 Neu-Ulm, Germany

MADE IN MALAYSIA