

Unigloves SOFT® NITRIL BLUE 300 (1410x*)

chemical tested · offline washed · blue
 long cuff · nitrile · single use gloves

Physical Specification:

Material: Nitrile (Flat)
 Colour: Blue
 Length: 300 mm
 Weight: 6.3 g
 Thickness: 0.12 mm
 Finish: Powder-free
 Fit: Ambidextrous
 Sterile: No

Order Codes | Sizing:

14102 | Small
 14103 | Medium
 14104 | Large
 14105 | Extra Large
 14106 | Extra Extra Large

Packaging:

Each case contains 10 boxes
 Each dispenser box contains 100 pairs

Classification:



For further
 information
 please click
 on the
 image



Features

- AQL 1.5
- Extended 300mm cuff
- Conforms to the hand
- Beaded cuff
- Textured Palms and Fingertips
- Touchscreen compatible

Benefits

- Meets medical standards for pinholes
- SOFT® NITRIL BLUE 300 gloves ensure no finger marks on metal or glass surfaces, maintaining pristine finishes
- Providing extra protection
- Designed to conform seamlessly to the hand, these gloves provide a snug and comfortable fit, enhancing dexterity and reducing hand fatigue
- The beaded cuff design ensures a secure fit and prevents the gloves from slipping off
- SOFT® NITRIL BLUE 300 gloves offer improved grip in wet or oily conditions, providing enhanced safety and control
- Users can operate mobile devices and tablets without the need to remove their gloves

Sectors

- Automotive maintenance and repair
- Electronics
- Glass
- Chemical manufacture and use
- Laboratory & pharmaceutical



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com



Issue No. 1

Prepared at Gillingham
 on 18/11/2025