

Unigloves CONTACT® (150X*)

offline washed · chemical tested · natural · medical grade
 medium weight · latex · examination gloves

Physical Specification:

Material: Latex (Flat)
 Colour: Natural
 Weight: 6.2 g
 Thickness: 0.11 mm
 Finish: Powder-free
 Fit: Ambidextrous
 Sterile: No

Order Codes | Sizing:

1501 | Extra Small
 1502 | Small
 1503 | Medium
 1504 | Large
 1505 | Extra Large

Packaging:

Each case contains 10 boxes
 Each dispenser box contains 100 pairs

Classification:



For further
 information
 please click
 on the
 image



Features

- Washed offline
- Tested to EN374 Type B
- AQL 1.5
- Free from powder
- Tested to EN1186
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Resistant to a range of chemicals
- Meets medical standards for pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves
- Extremely grippy and reliable on the outside, the glove feels soft against the skin on the inside and slides effortlessly onto the hand

Sectors

- Surgery
- Laboratories
- Clinics
- Chemical manufacture and use
- Dentistry
- Dental hygiene
- Emergency services




UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com



Issue No. 1

Prepared at Gillingham
 on 18/11/2025