

DERMA SKIN (170*) - THE SENSITIVE GLOVE

Disposable examination gloves made of latex

Latex • Powder free • Nature • Non-sterile •

Chemical protective glove EN 374 • Suitable for short-term food contact •

Protection against viruses, fungi and bacteria



Specifications:

Material: Latex
 Powder: Powder free
 Ingredients: Thiuram free
 Cuff: Beaded cuff
 Design: Ambidextrous
 Surface: Micro-roughened fingertips and palm
 Colour: Nature
 Quality control:
 EN 374 CAT III - Type B;
 EN ISO 21420;
 EN 455 Part 1-4;
 DIN EN ISO 13485:2021


Item number and size:

1701 | Extra Small
 1702 | Small
 1703 | Medium
 1704 | Large

Packaging:

100 pieces per box
 10 Boxen per carton

Classification:



Test parameter

Test criteria	Level	AQL
Dimensions	EN 455-2	Median value meets the requirements
Physical properties	EN 455-2	Median value meets the requirements
1000-ml-watertightness test	G-1	1.5
Visual inspection for major errors	G-1	2.5
Visual inspection for minor errors	G-1	4.0

Dimensions

REF	1701	1702	1703	1704
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9
Width	75 ± 5 mm	85 ± 5 mm	95 ± 5 mm	106 ± 5 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer (mm)
Cuff (at a distance of 40mm to the beaded edge)	0.11 ± 0.03
Palm (in the middle of the palm)	0.14 ± 0.03
Fingertips (at a distance of 15mm to the tip of the middle finger)	0.13 ± 0.03

Physical properties

Parameter	Before aging	After aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

CCQS Certification Services Ltd, Blanchardstown Corporate Park, D15 AKK1, Ireland
 (Notified Body 2834)

Test chemical	EN ISO 374-1:2018 + A1: 2018 Permeation level	EN ISO 374-4:2019 Degradation %
Sodium hydroxide 40% (K)	6	-55.11
Hydrogen peroxide 30% (P)	2	9.7
Formaldehyde 37% (T)	2	-12.4

Hint

The above information applies to the standard version.
 Subject to changes.



TOUCH SENSITIVITY



GRIP INTENSITY



SUPPLENESS &
SLIP-IN COMFORT



TEAR RESISTANCE



COMFORT & FIT



SKIN CARE



UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
 D-53842 Troisdorf
 ☎ +49 (0) 2241-9323-0
 ✉ unigloves@unigloves.de
 🌐 unigloves.de

