

# Vitality Soft Nitrile (GD005\*)

moisturising · offline washed · aqua blue · medical grade  
 light weight · nitrile · examination gloves

## Physical Specifications:

Material: Nitrile (latex free)  
 Colour: Aqua Blue  
 Length: 240 mm  
 Weight: 3.5 g  
 Thickness: 0.10 mm  
 Finish: Powder free, chlorinated  
 Fit: Ambidextrous  
 Sterile: No

## Order Codes and Size Information:

GD0051 | Extra small  
 GD0052 | Small  
 GD0053 | Medium  
 GD0054 | Large  
 GD0055 | Extra large

## Packaging:

Each dispenser box contains 100 gloves  
 Each case contains 10 dispenser boxes

## Classification:



For more about the **Vitality** range, please **click on the image**



## Features


- Washed offline
- Microfilm inner coating contains a combination of natural moisturisers
- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Touchscreen compatible

## Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Microfilm inner coating helps protect and nourish skin, and aids donning
- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

## Sectors

- Dentistry • Dental hygiene
- Floristry and gardening
- Hair and beauty
- Healthcare



**UNIGLOVES UK**  
 Registered in England No: 4010200  
 VAT registration 791817787

3 Ambley Green, Gillingham Business Park,  
 Gillingham, ME8 0NJ, UK  
 +44 (0)800 049 6602  
 enquiries@unigloves.co.uk  
 uk.unigloves.com

**UNIGLOVES GERMANY**  
 Arzt- & Klinikbedarf  
 Handelsgesellschaft mbH

Camp-Spich-Straße 71  
 D-53842 Troisdorf-Spich,  
 +49 (0) 2241-9323-0  
 unigloves@unigloves.de  
 unigloves.de

