

# Kooltouch (GM004\*)

offline washed · chemical tested · blue · medical grade  
medium weight · nitrile · examination gloves

**Physical Specifications:**

Material: Nitrile (latex free)  
Colour: Blue  
Length: 240 mm  
Weight: 4.7 g  
Thickness: 0.09 mm  
Finish: Powder free, chlorinated  
Fit: Ambidextrous  
Sterile: No

**Order Codes and Size Information:**

GM0041 | Extra small  
GM0042 | Small  
GM0043 | Medium  
GM0044 | Large  
GM0045 | Extra large

**Packaging:**

Each dispenser box contains 100 gloves  
Each case contains 10 dispenser boxes


**Classification:**

For more about the **Advanced** range, please **click on the image**

- Features**
- Washed offline
  - Tested to EN374
  - AQL 1.5 and EN455 parts 1-4
  - Tested to EN1186
  - Nitrile formulation conforms to the hand during use
  - Beaded cuff
  - Manufactured with textured palms and fingertips
  - Medium-weight construction
  - Touchscreen compatible

- Benefits**
- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
  - Resistant to a range of chemicals
  - Meets medical standards including pinholes
  - Suitable for contact with all food types
  - Nitrile formulation conforms to the hand during use for a more comfortable fit
  - Beaded cuff provides extra strength during donning and helps prevent liquid roll back
  - Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
  - Medium-weight construction provides a balance of strength and sensitivity
  - Suitable for use with mobile and tablet displays without removing gloves

- Sectors**
- Aerospace
  - Chemical manufacture and use
  - Dentistry • Dental hygiene
  - Electronics • Engineering
  - Food processing
  - Healthcare
  - Laboratory work
  - Manufacturing and assembly
  - Pharmaceutical work
  - Printing



**UNIGLOVES UK**  
Registered in England No: 4010200  
VAT registration 791817787

3 Ambley Green, Gillingham Business Park,  
Gillingham, ME8 0NJ, UK  
 +44 (0)800 049 6602  
 enquiries@unigloves.co.uk  
 uk.unigloves.com

**UNIGLOVES GERMANY**  
Arzt- & Klinikbedarf  
Handelsgesellschaft mbH

Camp-Spich-Straße 71  
D-53842 Troisdorf-Spich,  
 +49 (0) 2241-9323-0  
 unigloves@unigloves.de  
 unigloves.de

