

Pearl (GPO***)

12 colours · medical grade · medium weight nitrile · examination gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colours: Aqua Blue, Black, Blue, Burgundy, Green, Orange, Pink, Red, Sapphire, Violet, White, Yellow
 Length: 240 mm
 Weight: 3.8 g
 Thickness: 0.08 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GPO**1 | Extra small
 GPO**2 | Small
 GPO**3 | Medium
 GPO**4 | Large
 GPO**5 | Extra large

** 01 Blue / 02 White / 03 Black / 04 Green / 05 Pink / 06 Red / 07 Violet / 08 Sapphire / 09 Aqua Blue / 11 Yellow / 12 Burgundy / 13 Orange

Note: Not all colours come in Extra Small

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the Pearl range, please click on the image



Features

- Tested to EN374
- Available in twelve colours
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Ideal for colour-coding in food, cleaning and healthcare environments
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Care home work
- Chemical manufacture and use
- Cleaning and janitorial
- Dentistry • Dental hygiene • Emergency services
- Floristry and gardening • Food processing
- Hair and beauty • Healthcare
- Hospital and clinical duties
- Janitorial • Painting and decorating
- Piercing • Tattoo

Colour range:



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green, Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf-Spich,
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de

