

Unicare Blue Vinyl (GS009*)

blue · medium weight · vinyl · powdered
single use gloves

Physical Specifications:

Material: Vinyl (latex free)
Colour: Blue
Length: 240 mm
Weight: 4.5 g
Thickness: 0.06 mm
Finish: Lightly powdered
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GS0092 | Small
GS0093 | Medium
GS0094 | Large
GS0095 | Extra large

Packaging:

Each dispenser box contains 100 gloves
Each case contains 10 dispenser boxes

Classification:



For more about the Unicare range, please click on the image 

Features

- AQL 4.0
- Free from latex and accelerators
- Tested to EN1186
- Manufactured using high-quality vinyl
- Beaded cuff
- Medium-weight construction
- Lightly powdered

Benefits

- Meets industrial standards for pinholes
- Suitable for those with latex allergies (Type I sensitisation) and Type IV hypersensitivity
- Suitable for contact with non-fatty food types
- High-quality vinyl combines low cost with reduced risk of protein allergens
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Medium-weight construction provides a balance of strength and sensitivity
- USP grade corn starch powder improves ease of donning

Sectors

- Care home work
- Cleaning and janitorial
- Food processing
- General assembly
- Janitorial
- Local Authority work
- Painting and decorating
- Printing



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

📍 3 Ambley Green, Gillingham Business Park,
Gillingham, ME8 0NJ, UK
☎ +44 (0)800 049 6602
✉ enquiries@unigloves.co.uk
🌐 uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
D-53842 Troisdorf-Spich,
☎ +49 (0) 2241-9323-0
✉ unigloves@unigloves.de
🌐 unigloves.de

