

Unicare Yellow Vinyl (GS010*)

vinyl · powder free · lighter weight · non-sterile
single use examination gloves

Physical Specifications:

Material: Vinyl (latex free)
Colour: Yellow
Length: 240 mm
Weight: 4.5 g
Thickness: 0.06 mm
Finish: Powder free
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GS0102 | Small
GS0103 | Medium
GS0104 | Large
GS0105 | Extra large

Packaging:

Each dispenser box contains 100 gloves
Each case contains 10 dispenser boxes

Classification:



Food Safe
(not for use
with fatty foods)

Features

- Made from synthetic vinyl
- Unicare Yellow Vinyl gloves are light weight and dexterous
- AQL 4.0
- Latex free
- Accelerator free
- Powder free

Benefits

- Minimises the risk of skin issues
- Comfortable and easy to wear
- Industrial standard for pinholes
- Suitable for those with latex allergies (Type I) and delayed hypersensitivity (Type IV)
- Minimises particulate contamination

Sectors

- Food preparation
- Cleaning
- Care homes
- Food handling



For more about the Unicare range, please click on the image



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

📍 3 Ambley Green, Gillingham Business Park,
Gillingham, ME8 0NJ, UK
☎ +44 (0)800 049 6602
✉ enquiries@unigloves.co.uk
🌐 uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
D-53842 Troisdorf-Spich,
☎ +49 (0) 2241-9323-0
✉ unigloves@unigloves.de
🌐 unigloves.de

