



U

Unicare Red Vinyl (GS011*)

vinyl \cdot powder free \cdot lighter weight \cdot non-sterile single use examination gloves

Physical Specifications:

Material: Vinyl (latex free)

Colour: Red Length: 240 mm Weight: 4.5 g Thickness: 0.06 mm Finish: Powder free Fit: Ambidextrous Sterile: No

Order Codes and Size Information:

GS0112 | Small GS0113 | Medium GS0114 | Large GS0115 | Extra large

Packaging:

Each dispenser box contains 100 gloves Each case contains 10 dispenser boxes

Classification:







For more about the **Unicare** range, please click on the image



Features

- Made from synthetic vinyl
- Unicare Red Vinyl gloves are lighter weight and dexterous
- AQL 4.0
- Latex free
- Accelerator free
- Powder free

Benefits

- Minimises the risk of skin issues
- Comfortable and easy to wear
- Industrial standard for pinholes
- Suitable for those with latex allergies (Type I) and delayed hypersensitivity (Type IV)
- Minimises particulate contamination

♦ Sectors

- Food preparation
- Cleaning
- Care homes
- Food handling





UNIGLOVES UK

Registered in England No: 4010200 VAT registration 791817787

- 3 Ambley Green, Gillingham Business Park, Gillingham, ME8 ONJ, UK
- **2** +44 (0)800 049 6602
- enquiries@unigloves.co.uk
- uk.unigloves.com

UNIGLOVES GERMANY Arzt- & Klinikbedarf

Arzt- & Klinikbedarf Handelsgesellschaft mbH

- © Camp-Spich-Straße 71 D-53842 Troisdorf-Spich, 20 +49 (0) 2241-9323-0
- unigloves@unigloves.de

