

Unicare Flex Nitrile (GS021*)



blue · ultra light weight · nitrile
single use gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colour: Blue
 Length: 240 mm
 Weight: 2.2 g
 Thickness: 0.05 mm
 Finish: Powder free,
 chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GS0211 | Extra Small
 GS0212 | Small
 GS0213 | Medium
 GS0214 | Large
 GS0215 | Extra large

Packaging:

Each dispenser box
 contains 100 gloves
 Each case contains 10
 dispenser boxes

Classification:



For more about
 the Unicare
 range, please
 click
 on the
 image



Features


- AQL 4.0
- Free from powder
- Free from phthalates
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured fingertips
- Ultra-lightweight construction

Benefits

- Meets industrial standards for pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Compliant for contact with food
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured fingertips provide improved grip, particularly when working in wet or oily conditions
- Ultra-lightweight construction provides excellent tactile sensitivity with reduced cost

Sectors

- Care home work
- Cleaning and janitorial
- Food processing
- General assembly
- Local Authority work
- Painting and decorating
- Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green, Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf-Spich,
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de



Issue No. 7

Prepared at Gillingham
 on 22/11/22