

Unicare Black Nitrile (GS900*)

200 GLOVES

black · medical grade · light weight · nitrile examination gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colour: Black
 Length: 240 mm
 Weight: 3.5 g
 Thickness: 0.08 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GS9002 | Small
 GS9003 | Medium
 GS9004 | Large
 GS9005 | Extra large

Packaging:

Each dispenser box contains 200 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the Unicare range, please click on the image



Features


- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Free from phthalates
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured fingertips
- Lightweight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with food
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Lightweight construction provides good tactile sensitivity with reduced cost
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Care home work
- Cleaning and janitorial
- Dentistry • Dental hygiene
- Emergency services
- Food processing • General assembly
- Healthcare • Hospital and clinical duties
- Janitorial • Local Authority work
- Painting and decorating • Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de



Issue No. 1

Prepared at Gillingham
 on 25/07/25