

# Select Black Latex Extended Cuff (GT001\*)

extended cuff · offline washed · black · chemical tested  
medical grade · medium weight · latex examination gloves

**Physical Specifications:**

Material: Latex  
(offline chlorinated)  
Colour: Black  
Length: 300 mm  
Weight: 7.8 g  
Thickness: 0.13 mm  
Finish: Powder free, chlorinated  
Fit: Ambidextrous  
Sterile: No

**Order Codes and Size Information:**

GT0011 | Extra small  
GT0012 | Small  
GT0013 | Medium  
GT0014 | Large  
GT0015 | Extra large

**Packaging:**

Each dispenser box contains 100 gloves  
Each case contains 10 dispenser boxes

**Classification:**



For more about the **Select Black** range, please **click on the image**

**Features**

- Tested to EN374
- Extended cuff (300mm glove length)
- AQL 1.5 and EN455 parts 1-4
- Matt-black finish
- Free from powder
- Ambidextrous
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

**Benefits**

- Resistant to a range of chemicals
- Extended cuff helps prevent liquids getting into the glove when working in wet conditions and provides additional protection from chemical splashes
- Meets medical standards including pinholes
- Matt black finish masks blood and ink, and reduces hand and eye strain when working on complex designs
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

**Sectors**

- Hair and beauty
- Piercing
- Tattoo



**UNIGLOVES UK**  
Registered in England No: 4010200  
VAT registration 791817787

📍 3 Ambley Green, Gillingham Business Park, Gillingham, ME8 0NJ, UK  
☎ +44 (0)800 049 6602  
✉ enquiries@unigloves.co.uk  
🌐 uk.unigloves.com

**UNIGLOVES GERMANY**  
Arzt- & Klinikbedarf  
Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71, D-53842 Troisdorf-Spich,  
☎ +49 (0) 2241-9323-0  
✉ unigloves@unigloves.de  
🌐 unigloves.de

