


Select Black Nitrile (GT003*)

black · chemical tested · medical grade · medium weight nitrile · examination gloves



Physical Specifications:

Material: Nitrile (offline chlorinated)
 Colour: Black
 Length: 240 mm
 Weight: 3.8 g
 Thickness: 0.08 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GT0031 | Extra small
 GT0032 | Small
 GT0033 | Medium
 GT0034 | Large
 GT0035 | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the **Select Black** range, please **click on the image**

Features


- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Manufactured with an enhanced fingertips grip
- Matt-black finish
- Free from latex
- Free from powder
- Ambidextrous
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Medium-weight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Grippy fingertips means less effort to hold instruments, particularly when working in wet or oily conditions
- Matt black finish masks blood and ink, and reduces hand and eye strain when working on complex designs
- Suitable for those with latex allergies (Type I sensitisation)
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Hair and beauty
- Piercing
- Tattoo



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

📍 3 Ambley Green, Gillingham Business Park, Gillingham, ME8 0NJ, UK
 ☎ +44 (0)800 049 6602
 ✉ enquiries@unigloves.co.uk
 🌐 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71, D-53842 Troisdorf-Spich,
 ☎ +49 (0) 2241-9323-0
 ✉ unigloves@unigloves.de
 🌐 unigloves.de

