



SELECT Black 300 (GT2430X*)

offline washed · black · chemical tested · medical grade
medium weight · latex examination gloves

Physical Specification:

Material: Latex (Flat)
Colour: Black
Weight: 7.8 g
Thickness: 0.13 mm
Finishing: Powder-free
Fit: Ambidextrous
Sterile: No

Order Codes | Sizing:

GT24301 | Extra Small
GT24302 | Small
GT24303 | Medium
GT24304 | Large
GT24305 | Extra Large

Packaging:

Each case contains 10 boxes
Each dispenser box
contains 100 pairs

Classification:



For further
information
please click
on the
image



Features

- Tested to EN374 Type B
- Extended cuff (300mm glove length)
- AQL 1.5
- Matt-black finish
- Free from powder
- Ambidextrous
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Extended cuff helps prevent liquids getting into the glove when working in wet conditions and provides additional protection from chemical splashes
- Meets medical standards for pinholes
- Matt black finish masks blood and ink, and reduces hand and eye strain when working on complex designs
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Hair and beauty
- Piercing
- Tattoo



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
de.unigloves.com

