

# PRO.TECT Orange Diamond HD (GA005\*)

silicone free · diamond textured · chemical tested  
orange · ultra heavy weight · nitrile · single use gloves

## Physical Specifications:

Material: Nitrile (latex free)  
Colour: Orange  
Length: 240 mm  
Weight: 7.6 g  
Thickness:  
0.13 mm (flat area)  
0.22 mm (raised area)  
Finish: Powder free,  
chlorinated  
Fit: Ambidextrous  
Sterile: No

## Order Codes and Size Information:

GA0052 | Small  
GA0053 | Medium  
GA0054 | Large  
GA0055 | Extra large

## Packaging:

Each dispenser box contains 100 gloves  
Each case contains 10 dispenser boxes

## Classification:



For more about the PRO.TECT range, please click on the image 



## Features


- Free from silicone
- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Ambidextrous
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with diamond textured palms and fingers (and back of hand)
- Ultra-heavyweight construction
- Touchscreen compatible

## Benefits

- Silicone free formulation means no finger marks on metal and glass surfaces
- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Diamond textured palms and fingers provides improved grip, particularly when working in wet or oily conditions
- Ultra-heavyweight construction is very tough, durable and long lasting
- Suitable for use with mobile and tablet displays without removing gloves

## Sectors

- Automotive • Automotive maintenance and repair
- Chemical manufacture and use • Construction
- Electronics • Engineering • Floristry and gardening
- General assembly • Repair • Glass • Inspection
- Laboratory work • Light duty maintenance • Utility
- Maintenance • Manufacturing and assembly
- Painting and decorating



full palm & finger raised grip pattern on both sides of glove



**UNIGLOVES UK**  
Registered in England No: 4010200  
VAT registration 791817787

 3 Ambley Green, Gillingham Business Park,  
Gillingham, ME8 0NJ, UK  
 +44 (0)800 049 6602  
 enquiries@unigloves.co.uk  
 uk.unigloves.com

**UNIGLOVES GERMANY**  
Arzt- & Klinikbedarf  
Handelsgesellschaft mbH

 Camp-Spich-Straße 71  
D-53842 Troisdorf-Spich,  
 +49 (0) 2241-9323-0  
 unigloves@unigloves.de  
 unigloves.de

