

Kooltouch (GM004*)

offline washed · chemical tested · blue · medical grade
medium weight · nitrile · examination gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colour: Blue
 Length: 240 mm
 Weight: 4.4 g
 Thickness: 0.09 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GM0041 | Extra small
 GM0042 | Small
 GM0043 | Medium
 GM0044 | Large
 GM0045 | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the **Advanced** range, please **click on the image**



Features


- Washed offline
- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace
- Chemical manufacture and use
- Dentistry • Dental hygiene
- Electronics • Engineering
- Food processing
- Healthcare
- Laboratory work
- Manufacturing and assembly
- Pharmaceutical work
- Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com

