

Stronghold (GM006*)

AQL 0.65 · chemical tested · chemotherapy drug tested · purple medical grade · heavy weight · nitrile · examination gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colour: Purple
 Length: 240 mm
 Weight: 6.2 g
 Thickness: 0.12 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GM0061 | Extra small
 GM0062 | Small
 GM0063 | Medium
 GM0064 | Large
 GM0065 | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the **Advanced** range, please click on the image



Features


- Free from silicone
- Tested to EN374
- Tested to ASTM D6978-05
- AQL 0.65 and EN455 parts 1-4
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Heavy-weight construction
- Touchscreen compatible

Benefits

- Silicone free formulation means no finger marks on metal and glass surfaces
- Resistant to a range of chemicals
- Resistant to a range of cytotoxic drugs when tested to ASTM D6978-05
- Meets medical standards + ultra-low hole level (exceeds medical standards for pinholes)
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Heavy-weight construction is tough, durable and long lasting
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace
- Chemical manufacture and use
- Dentistry • Dental hygiene
- Electronics • Engineering
- Food processing
- Healthcare • Laboratory work
- Manufacturing and assembly
- Pharmaceutical work • Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com

