

Stronghold+ (GM007*)

extended cuff · AQL 0.65 · enhanced chemical splash protection
 silicone free · purple · medical grade · heavy weight · nitrile
 unique Unigloves formulation · single use gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colour: Purple
 Length: 300 mm
 Weight: 7.3 g
 Thickness: 0.12 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GM0071 | Extra small
 GM0072 | Small
 GM0073 | Medium
 GM0074 | Large
 GM0075 | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the **Advanced** range, please click on the image 



Features


- Free from silicone
- Tested to EN374
- Tested to ASTM D6978-05
- Extended cuff (300mm glove length)
- AQL 0.65 and EN455 parts 1-4
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingers
- Heavy-weight construction
- Touchscreen compatible

Benefits

- Silicone free formulation means no finger marks on metal and glass surfaces
- Resistant to a range of chemicals
- Resistant to a range of cytotoxic drugs when tested to ASTM D6978-05
- Extended cuff helps prevent liquids getting into the glove when working in wet conditions and provides additional protection from chemical splashes
- Meets medical standards + ultra-low hole level (exceeds medical standards for pinholes)
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingers provide improved grip, particularly when working in wet or oily conditions
- Heavy-weight construction is tough, durable and long lasting
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace • Chemical manufacture and use • Dentistry
- Dental hygiene • Electronics • Engineering • Food processing
- Healthcare • Laboratory work • Manufacturing and assembly
- Pharmaceutical work • Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com

