

Product specification sheet

Revision dated October 2025

1.	Product code	GS001*									
2.	Product description	Unicare Safetec Latex									
3.	Material	Latex									
4.	Inspection parameters	Criteria		Inspection level		AQL		Standard			
		Dimensions		S-2		4.0		EN455-2			
		Physical properties		S-2		4.0		EN455-2			
		Water tight test (1,000ml)		G-1		1.5		EN455-1			
		Visual defects		G-1		2.5		In-house			
5.	Glove dimensions	Dimensions		Size (mm)							
				XS		S	M		L	XL	
		Product code		GS0011	GS0012	GS0013		GS0014		GS0015	
		Palm width		≤80	80 ± 10	95 ± 10		110 ± 10		≥ 110	
		Length		≥ 240	≥ 240	≥ 240		≥ 240		≥ 240	
6.	Glove thickness	Point of thickness measurement						Nominal single wall (mm)			
		Cuff (at beaded end)						0.08			
		Palm (at centre of palm)						0.10			
		Fingertip (at fingertip)						0.12			
7.	Physical properties	Parameters		Before aging				After aging			
		Median force at break (N)		Min 6.0				N/A			
		Tensile strength (MPa)		Min 18.0				Min 14.0			
		Ultimate elongation (%)		Min 650				Min 500			
8.	Proteins	< 50 µg/g extractable latex proteins									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Microroughened									
13.	Colour	Natural									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	5.7g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN ISO 21420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.