

# Unicare Blue Nitrile (GS004\*)

**200 GLOVES**

blue · medical grade · light weight · nitrile examination gloves

## Physical Specifications:

Material: Nitrile (latex free)  
 Colour: Blue  
 Length: 240 mm  
 Weight: 3.5 g  
 Thickness: 0.08 mm  
 Finish: Powder free, chlorinated  
 Fit: Ambidextrous  
 Sterile: No

## Order Codes and Size Information:

GS0041 | Extra small  
 GS0042 | Small  
 GS0043 | Medium  
 GS0044 | Large  
 GS0045 | Extra large

## Packaging:

Each dispenser box contains 200 gloves  
 Each case contains 10 dispenser boxes

## Classification:



For more about the Unicare range, please click on the image



## Features


- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Free from phthalates
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured fingertips
- Lightweight construction
- Touchscreen compatible

## Benefits

- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with food
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Lightweight construction provides good tactile sensitivity with reduced cost
- Suitable for use with mobile and tablet displays without removing gloves

## Sectors

- Care home work
- Cleaning and janitorial
- Dentistry • Dental hygiene
- Emergency services
- Food processing • General assembly
- Healthcare • Hospital and clinical duties
- Janitorial • Local Authority work
- Painting and decorating • Printing




**UNIGLOVES UK**  
 Registered in England No: 4010200  
 VAT registration 791817787

3 Ambley Green,  
 Gillingham Business Park,  
 Gillingham, ME8 0NJ, UK  
 +44 (0)800 049 6602  
 enquiries@unigloves.co.uk  
 uk.unigloves.com

**UNIGLOVES GERMANY**  
 Arzt- & Klinikbedarf  
 Handelsgesellschaft mbH

Camp-Spich-Straße 71  
 D-53842 Troisdorf,  
 Germany  
 +49 (0) 2241-9323-0  
 unigloves@unigloves.de  
 unigloves.de



Issue No. 6

Prepared at Gillingham  
 on 09/06/25