

Unicare Clear Vinyl (GS006*-A)

clear · medical grade · medium weight · vinyl examination gloves

Physical Specifications:

Material: Vinyl (latex free)
 Colour: Clear
 Length: 240 mm
 Weight: 4.5 g
 Thickness: 0.06 mm
 Finish: Powder free
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GS0062-A | Small
 GS0063-A | Medium
 GS0064-A | Large
 GS0065-A | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the Unicare range, please click on the image



Features

- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Free from latex and accelerators
- Free from powder
- Tested to EN1186
- Ambidextrous
- Manufactured using high-quality vinyl
- Beaded cuff
- Medium-weight construction

Benefits

- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Suitable for those with latex allergies (Type I sensitisation) and Type IV hypersensitivity
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with non-fatty food types
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- High-quality vinyl combines low cost with reduced risk of protein allergens
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Medium-weight construction provides a balance of strength and sensitivity

Sectors

- Care home work
- Cleaning and janitorial
- Dentistry
- Dental hygiene
- Emergency services
- Food processing
- General assembly
- Healthcare
- Hospital and clinical duties
- Janitorial
- Local Authority work
- Painting and decorating
- Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de

