

Unicare Blue TPE (GS023*)

blue · light weight · TPE
single use gloves

Physical Specifications:

Material: TPE
Colour: Blue
Length: 255 mm
Weight: 2.0 g
Thickness: -
Finish: Powder free,
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GS0232 | Small
GS0233 | Medium
GS0234 | Large
GS0235 | Extra large

Packaging:

Each dispenser box
contains 200 gloves
Each case contains 10
dispenser boxes

Classification:



For more about
the **Unicare**
range, please
click
on the
image



Features

- AQL 4.0
- Made from TPE (Thermoplastic Elastomer)
- Free from latex and accelerators
- Free from powder
- Stretch-to-hand fit
- Free from phthalates
- Tested to EN1186
- Lightweight construction
- Touchscreen compatible

Benefits

- Meets industrial standards for pinholes
- TPE is a flexible, comfortable and affordable material - the ideal alternative to vinyl gloves
- Suitable for those with latex allergies (Type I sensitisation) and delayed hypersensitivity (type IV)
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Stretch-to-hand fit for precise, dextrous movements
- Suitable for contact with food
- Suitable for contact with all food types
- Lightweight construction provides good tactile sensitivity with reduced cost
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Cleaning industry and household
- Food industry
- Janitorial
- General purpose
- Painting & DIY



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
unigloves.de

