

Product specification sheet

Revision dated October 2025

1.	Product code	GP012*									
2.	Product description	Burgundy Pearl Nitrile									
3.	Material	Nitrile									
4.	Inspection parameters	Criteria		Inspection level		AQL		Standard			
		Dimensions		S-2		4.0		EN455-2			
		Physical properties		S-2		4.0		EN455-2			
		Water tight test (1,000ml)		G-1		1.5		EN455-1			
		Visual defects		G-1		2.5		In-house			
5.	Glove dimensions	Dimensions		Size (mm)							
				XS		S	M	L	XL		
		Product code		GP0121	GP0122	GP0123	GP0124	GP0125			
		Palm width		≤80	80 ± 10	95 ± 10	110 ± 10	≥ 110			
		Length		≥ 240	≥ 240	≥ 240	≥ 240	≥ 240			
6.	Glove thickness	Point of thickness measurement					Nominal single wall (mm)				
		Cuff (at beaded end)					0.05				
		Palm (at centre of palm)					0.07				
		Fingertip (at fingertip)					0.10				
7.	Physical properties	Parameters		Before aging			After aging				
		Median force at break (N)		Min 6.0			Min 6.0				
		Tensile strength (MPa)		Min 14.0			Min 14.0				
		Ultimate elongation (%)		Min 500			Min 400				
8.	Proteins	None									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Fingertip textured									
13.	Colour	Burgundy									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	3.5g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN ISO 21420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.