

Product specification sheet

Revision dated October 2025

1.	Product code	GS002*									
2.	Product description	Unicare Classic Latex powdered									
3.	Material	Latex									
4.	Inspection parameters	Criteria		Inspection level		AQL		Standard			
		Dimensions		S-2		4.0		EN455-2			
		Physical properties		S-2		4.0		EN455-2			
		Water tight test (1,000ml)		G-1		1.5		EN455-1			
		Visual defects		G-1		2.5		In-house			
5.	Glove dimensions	Dimensions		Size (mm)							
				XS		S	M	L	XL		
		Product code		GS0021	GS0022	GS0023	GS0024	GS0025			
		Palm width		≤80	80 ± 10	95 ± 10	110 ± 10	≥ 110			
		Length		≥ 240	≥ 240	≥ 240	≥ 240	≥ 240			
6.	Glove thickness	Point of thickness measurement					Nominal single wall (mm)				
		Cuff (at beaded end)					0.08				
		Palm (at centre of palm)					0.09				
		Fingertip (at fingertip)					0.10				
7.	Physical properties	Parameters		Before aging			After aging				
		Median force at break (N)		Min 6.0			N/A				
		Tensile strength (MPa)		Min 18.0			Min 14.0				
		Ultimate elongation (%)		Min 650			Min 500				
8.	Proteins	> 50 µg/g extractable latex proteins									
9.	Powder residue	Maximum 10.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Microroughened									
13.	Colour	Natural									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	5.0g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN ISO 21420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓					

The above data represents the product norm, which may be changed without notice and is intended as a guide only.