

BioTouch (GM008* & GM009*)

biodegradable · AQL 1.5 · chemical tested · 2 colours
 medical grade · medium weight · nitrile · examination gloves



Physical Specifications:

Material: Nitrile (latex free)
 Colours: Violet blue, Black
 Length: 240 mm
 Weight: 3.5 g
 Thickness: 0.07 mm
 Finish: Powder free
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

Violet blue

GM0081 | Extra Small
 GM0082 | Small
 GM0083 | Medium
 GM0084 | Large
 GM0085 | Extra large

Black

GM0091 | Extra Small
 GM0092 | Small
 GM0093 | Medium
 GM0094 | Large
 GM0095 | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about the **Advanced** range, please **click on the image**



Features


- >81% biodegradation within 490 days*
- Contains innovative biodegradable technology
- Tested to EN374 Type B
- Tested to ASTM D5511 and ASTM D5526
- AQL 1.5 and EN455 parts 1-4
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured fingertips
- Touchscreen compatible

Benefits

- Environmentally friendly, much faster organic break down compared to regular nitrile, which can take over 100 years
- The technology means the biodegradation process only starts when disposed of in an active landfill environment
- Resistant to a range of chemicals
- Tested to leading biodegradation standards
- Meets medical standards including pinholes
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Food
- Automotive
- Aerospace
- Chemical manufacture and use
- Dentistry
- Dental hygiene
- Electronics
- Engineering
- Food processing
- Healthcare
- Laboratory work
- Manufacturing and assembly
- Pharmaceutical work
- Printing
- Tattoo



* Above results are based on a real-time study. The actual biodegradation rates will vary depending on the landfill conditions and the biological activity of microorganisms surround the gloves.

UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com

