

Unicare Classic Latex (GS002*)

natural · medium weight · latex · powdered
single use gloves

Physical Specifications:

Material: Latex
 Colour: Natural
 Length: 240 mm
 Weight: 5.0 g
 Thickness: 0.09 mm
 Finish: Lightly powdered,
 chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GS0021 | Extra small
 GS0022 | Small
 GS0023 | Medium
 GS0024 | Large
 GS0025 | Extra large

Packaging:

Each dispenser box
 contains 100 gloves
 Each case contains 10
 dispenser boxes

Classification:



For more about
 the Unicare
 range, please
 click
 on the
 image



Features


- AQL 1.5 and EN455 parts 1-4
- Beaded cuff
- Manufactured with textured palms and fingertips
- Lightweight construction
- Lightly powdered

Benefits

- Meets medical standards including pinholes
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Lightweight construction provides good tactile sensitivity with reduced cost
- USP grade corn starch powder improves ease of donning

Sectors

- Care home work
- Cleaning and janitorial
- General assembly
- Janitorial
- Local Authority work
- Painting and decorating
- Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de

