

Unigloves COMFORT® 300 (60X*)

extended cuff · natural · offline washed · chemical tested
 medical grade · medium weight · latex examination gloves

Physical Specification:

Material: Latex (Flat)
 Sterile: No
 Weight: 7.8 g
 Powder: Powder-free
 Colour: Natural
 Thickness: 0.11 mm
 Fit: Ambidextrous

Order Codes | Sizing:

602 | Small
 603 | Medium
 604 | Large

Packaging:

Each case contains 10 boxes
 Each dispenser box contains 100 pairs

Classification:



For further information
 please click
 on the image



Features

- Washed offline
- Tested to EN374 Type B
- AQL 1.5
- Free from powder
- Tested to EN1186
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Extended cuff helps prevent liquids getting into the glove when working in wet conditions and provides additional protection from chemical splashes
- Resistant to a range of chemicals
- Meets medical standards for pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves"

Sectors

- Aerospace • Chemical manufacture and use
- Dentistry • Dental hygiene • Electronics
- Engineering • Food processing • Healthcare
- Laboratory work • Manufacturing and assembly
- Pharmaceutical work • Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com



Issue No. 1

Prepared at Gillingham
 on 18/11/2025