

Unigloves Soft Nitrile Blue (80X0*)

Chemical tested · blue
medium weight nitrile single use gloves

Physical Specification:

Material: Nitrile (Flat)
Colour: Blue
Weight: 4.5 g
Thickness: 0.08 mm
Finish: Powder-free
Fit: Ambidextrous
Sterile: No

Order Codes | Sizing:

8020 | Small
8030 | Medium
8040 | Large
8050 | Extra Large

Packaging:

Each case contains 10 boxes
Each dispenser box contains 100 pairs

Classification:



For further
information
please click
on the
image



Features

- Tested to EN374
- AQL 1.5
- Free from powder
- Free from phthalates
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Lightweight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Meets medical standards for pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with food
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Lightweight construction provides good tactile sensitivity with reduced cost
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Care home work • Cleaning and janitorial
- Dentistry • Dental hygiene • Emergency services
- Food processing • General assembly
- Healthcare • Hospital and clinical duties
- Janitorial • Local Authority work
- Painting and decorating • Printing



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
de.unigloves.com



Issue No. 1

Prepared at Gillingham
on 18/11/2025