

Pearl Violet (GP007*)

violet · medical grade · medium weight
nitrile · examination gloves

Physical Specifications:

Material: Nitrile (latex free)
 Colour: Violet
 Length: 240 mm
 Weight: 3.5 g
 Thickness: 0.07 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GP0071 | Extra small
 GP0072 | Small
 GP0073 | Medium
 GP0074 | Large
 GP0075 | Extra small

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



For more about our **Nitrile** range, please **click on the image** 



Features


- Tested to EN374
- Available in sixteen colours
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Ideal for colour-coding in food, cleaning and healthcare environments
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Care home work
- Chemical manufacture and use
- Cleaning and janitorial
- Dentistry • Dental hygiene
- Emergency services
- Floristry and gardening
- Food processing • Hair and beauty
- Healthcare • Hospital and clinical duties
- Janitorial • Painting and decorating
- Piercing • Tattoo



UNIGLOVES USA

UG Healthcare USA Inc.
 Costa Mesa, CA 92626

(800) 328-7981
 sales@ughealthcare.com
 us.unigloves.com



Issue No. 1

Prepared at Gillingham
 on 23/09/26