

Cobalt Pearl (GP101*) Einweg-Untersuchungshandschuhe aus Nitril Nitril • Puderfrei • Cobalt • Lebensmittelecht • Unsteril Chemikalienschutzhandschuh EN 374 Schutz vor Viren, Pilzen und Bakterien

Spezifikationen:

Material: Nitril
 Puder: Ungepudert
 Inhaltsstoffe: Thiuramfrei,
 keine Proteine enthalten
 Manschette: Umbördelt
 Design: Beidhändig
 Oberfläche: Angeraute
 Fingerspitzen
 Farbe: Cobalt
 Qualitätssicherung:
 EN 374 CAT III - Typ C;
 EN ISO 21420;
 EN 455 Teil 1-4;
 DIN EN ISO 13485:2021

Artikelnummer und Größe:

GP1011 | Extra Small
 GP1012 | Small
 GP1013 | Medium
 GP1014 | Large
 GP1015 | Extra large

Verpackung:

100 Stück pro Box
 10 Boxen pro Karton

Klassifikation:



Prüfparameter

Prüfkriterien	Prüfstufe	AQL
Maße	EN 455-2 Medianwert entspricht den Anforderungen	
Physikalische Eigenschaften	EN 455-2 Medianwert entspricht den Anforderungen	
1000-ml-Wasserdichtheitstest	G-1	1,5
Sichtprüfung auf größere Fehler	G-1	2,5
Sichtprüfung auf kleinere Fehler	G-1	4,0

Maße

Art.-Nr.	GP1011	GP1012	GP1013	GP1014	GP1015
Größe	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Breite	70 ± 10 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	120 ± 10 mm
Länge	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Materialstärke

Messpunkte	Eine Lage
Handfläche (in der Mitte der Handfläche)	0,09 ± 0,03 mm

Physikalische Eigenschaften

Parameter	Vor der Alterung	Nach der Alterung
Reißkraft im Median (N)	min. 6,0	min. 6,0

Chemikalienschutz

Das Produkt entspricht der PSA-Verordnung (EU) 2016/425

Baumusterprüfung (Modul B) & Modul C2
 CCQS Certification Services Ltd, Blanchardstown Corporate Park, D15 AKK1, Ireland
 (Notified Body 2834)

Prüfchemikalie	EN ISO 374-1:2016 Permeationslevel	EN374-4:2019 Degradation %
Natriumhydroxid 40% (K)	6	-5,3

Hinweis

Die vorstehenden Angaben gelten für die Standardausführung.
 Änderungen vorbehalten.




UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com





Cobalt Pearl (GP101*)

Disposable examination gloves made of nitrile

Nitrile • Powder-free • Cobalt • Suitable for short-term food contact

Non-sterile • Chemical protective glove EN 374

Protection against viruses, fungi and bacteria

Specifications:

Material: Nitrile
 Powder: Powder-free
 Ingredients: Thiuram free, no proteins contained
 Cuff: Beaded cuff
 Design: Ambidextrous
 Surface: Finger textured
 Colour: Cobalt
 Quality assurance:
 EN 374 CAT III - Type C;
 EN ISO 21420;
 EN 455 Part 1-4;
 DIN EN ISO 13485:2021

Item number and size:

GP1011 | Extra Small
 GP1012 | Small
 GP1013 | Medium
 GP1014 | Large
 GP1015 | Extra large

Packaging:

100 pieces per box
 10 boxes per carton

Classification:



Test parameter

Test criteria	Level	AQL
Dimensions	EN 455-2	Median value meets the requirements
Physical properties	EN 455-2	Median value meets the requirements
1000-ml-Watertightness test	G-1	1.5
Visual inspection for major defects	G-1	2.5
Visual inspection for minor defects	G-1	4.0

Dimensions

Item no.	GP1011	GP1012	GP1013	GP1014	GP1015
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	70 ± 10 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	120 ± 10 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer
Palm (in the middle of the palm)	0.09 ± 0.03 mm

Physical properties

Parameter	Before aging	After aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

CCQS Certification Services Ltd, Blanchardstown Corporate Park, D15 AKK1, Ireland
 (Notified Body 2834)

Test chemical	EN ISO 374-1:2016 Permeation level	EN374-4:2019 Degradation %
Sodium hydroxide 40% (K)	6	-5.3

Hint

The above information applies to the standard version.
 Subject to changes.



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com

