

Unicare Soft Blue Nitrile (GS004*)

200 GLOVES

blue · medical grade · light weight · nitrile examination gloves

Physical Specifications:

Material: Nitrile (latex free)
Colour: Blue
Length: 240 mm
Weight: 3.5 g
Thickness: 0.08 mm
Finish: Powder free, chlorinated
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GS0041 | Extra small
GS0042 | Small
GS0043 | Medium
GS0044 | Large
GS0045 | Extra large

Packaging:

Each dispenser box contains 200 gloves
Each case contains 10 dispenser boxes

Classification:



For more about our **Nitrile** range, please **click on the image**



Features


- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Free from phthalates
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured fingertips
- Lightweight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with food
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Lightweight construction provides good tactile sensitivity with reduced cost
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Care home work
- Cleaning and janitorial
- Dentistry • Dental hygiene
- Emergency services
- Food processing • General assembly
- Healthcare • Hospital and clinical duties
- Janitorial • Local Authority work
- Painting and decorating • Printing



UNIGLOVES USA

UG Healthcare USA Inc.
Costa Mesa, CA 92626

(800) 328-7981
sales@ughealthcare.com
us.unigloves.com



Issue No. 1

Prepared at Gillingham
on 03/12/25