

UNICARE BLACK LATEX (GS290*)

Einweg-Untersuchungshandschuhe aus Latex

Latex • puderfrei • schwarz • unsteril • Chemikalienschutzhandschuh EN 374
Schutz vor Viren, Pilzen und Bakterien

Spezifikationen:

Material: Hochwertiges Naturlatex
Puder: Ungepudert
Manschette: Umbördelt
Design: Beidhändig
Oberfläche: Mikrogeraute Handfläche
Farbe: Schwarz
Qualitätssicherung:
EN 374 CAT III - Typ B,
EN 420; EN 455 Teil 1-4;
EN ISO 13485:2016


Artikelnummer und Größe:

GS2902 | Small
GS2903 | Medium
GS2904 | Large
GS2905 | Extra Large

Verpackung:

100 Stück pro Box
10 Boxen pro Karton

Klassifikation:



Für weitere Informationen klicken Sie bitte auf das Bild



Prüfparameter

Prüfkriterien	Prüfstufe	AQL
Maße	S-2	4,0
Physikalische Eigenschaften	S-2	4,0
1000-ml-Wasserdichtheitstest	G-1	1,5
Sichtprüfung auf Fehler	G-1	2,5

Maße

Art.-Nr.	GS2902	GS2903	GS2904	GS2905
Größe	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Breite	84 ± 3 mm	94 ± 3 mm	105 ± 3 mm	113 ± 3 mm
Länge	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Materialstärke

Meßpunkte	Eine Lage (mm)
Handfläche (in der Mitte der Handfläche)	0.11 ± 0.02
Fingerspitze (im Abstand von 15mm zur Spitze des Mittelfingers)	0.12 ± 0.02

Physikalische Eigenschaften

Parameter	Vor der Alterung	Nach der Alterung
Zugfestigkeit (N)	min. 6.0	min. 6.0

Chemikalienschutz

Das Produkt entspricht der PSA-Verordnung (EU) 2016/425

Type-Examination (Modul B) & Modul C2

SATRA Technology Europe, Bracetown Business Park, Clonee, D15 YN2, Ireland.

(Notified Body 2777)

Prüfchemikalie	EN ISO 374-1:2016+A1: 2018/Typ B Permeationslevel	EN ISO 374-4:2019 Degradation %
40% Sodium Hydroxide (K)	6	-76.0
96% Sulphuric Acid (L)	3	72.1
65% Nitric Acid (M)	2	21.8
30% Hydrogen Peroxide (P)	6	-8.6
37% Formaldehyde (T)	3	-40.8

Hinweis

Die vorstehenden Angaben gelten für die Standardausführung. Änderungen vorbehalten.



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
unigloves.de



UNICARE BLACK LATEX (GS290*)

Disposable Latex Examination Gloves

Latex • powder free • black • non-sterile • Chemical protective glove EN 374 protection against viruses, fungi and bacteria

Specifications:

Material: Natural High Grade Rubber Latex
 Powder: powder free
 Cuff: beaded cuff
 Design: Ambidextrous
 Surface: Micro-roughened palm
 Colour: black
 Quality assurance:
 EN 374 CAT III - Type B,
 EN 420; EN 455 Part 1-4;
 EN ISO 13485:2016


Item No. and Size:

GS2902 | Small
 GS2903 | Medium
 GS2904 | Large
 GS2905 | Extra Large

Packaging:

100 Pieces per Box
 10 Boxes per Carton

Classification:



For more Informations
 Klick on the Image



Check parameters

Characteristics	Level	AQL
Dimensions	S-2	4,0
Physical Properties	S-2	4,0
1000-ml-Watertightness test	G-1	1,5
Visual inspection for errors	G-1	2,5

Dimensions

Item No.	GS2902	GS2903	GS2904	GS2905
Size	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	84 ± 3 mm	94 ± 3 mm	105 ± 3 mm	113 ± 3 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Thickness

Measuring Points	One Layer (mm)
Palm (in the middle of the palm)	0.11 ± 0.02
Fingertips (distance 15 mm to top of middle finger)	0.12 ± 0.02

Physical Properties

Parameter	Before Aging	After Aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type-Examination (Modul B) & Modul C2

SATRA Technology Europe, Bracetown Business Park, Clonee, D15 YN2P. Ireland
 (Notified Body 2777)

Test chemical	EN ISO 374-1:2016+A1: 2018/Typ B Permeation level	EN ISO 374-4:2019 Degradation %
40% Sodium Hydroxide (K)	6	-76.0
96% Sulphuric Acid (L)	3	72.1
65% Nitric Acid (M)	2	21.8
30% Hydrogen Peroxide (P)	6	-8.6
37% Formaldehyde (T)	3	-40.8

Hint

The above information apply to the standard version.
 Subject to changes.



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de

