

UNICARE WHITE NITRILE (GS700*)

Disposable examination gloves made of nitrile

Nitrile • Latex free • Powder-free • White • Non-sterile • Single chlorinated
Chemical protective glove EN 374 • Suitable for short-term food contact
Protection against viruses, fungi and bacteria

Specifications:

Material: Nitrile
Powder: Powder-free
Ingredients: Thiuram free, no proteins contained
Cuff: Beaded cuff
Design: Ambidextrous
Surface: Finger textured
Colour: White
Quality assurance:
EN 374 CAT III - Type C;
EN ISO 21420;
EN 455 Part 1-4;
DIN EN ISO 13485:2021

Item Number and Size:

GS7002 | Small
GS7003 | Medium
GS7004 | Large
GS7005 | Extra Large

Packaging:

200 pieces per box
10 boxes per carton

Classification:



For more information click on the Image 



Test parameter

Test criteria	Level	AQL
Dimensions	EN 455-2 Median value meets the requirements	
Physical properties	EN 455-2 Median value meets the requirements	
1000-ml-Watertightness test	G-1	1.5
Visual inspection for major defects	G-1	2.5
Visual inspection for minor defects	G-1	4.0

Dimensions

Item no.	GS7002	GS7003	GS7004	GS7005
Size	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	min. 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer
Palm (in the middle of the palm)	0.08 ± 0.03 mm

Physical Properties

Parameter	Before aging	After aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

CCQS Certification Services Limited, Block 1 Blanchardstown Corporate Park,
Ballycoolin Road, Blanchardstown, D15 AKK1, Ireland
(Notified Body 2834)

Test chemical	EN ISO 374-1:2016+A1:2018/Typ C Permeation level	EN ISO 374-4:2019 Degradation %
Sodium hydroxide 40% (K) 6		-5.3

Hint

The above information applies to the standard version.
Subject to changes.



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
unigloves.de

